

Specifications

Pocket sized and feather-weight for portable use in a large number of applications. It is utilize user-replaceable electrodes

Measuring range: 0.00-14.00pH

Accuracy: ± 0.1 PH

Resolution: 0.01PH

Operating temperature : 0C- 50C

Automatic temperature compensation: 0C- 50C

Plug: BNC+1m cable

Dimensions: 120mm x 26mm x 26mm

Weight: 70g

Operation

- 1.Remove the protective cap of electrode.
- 2.First rinse the electrode with distilled water, and suck it with filter-paper.
- 3.Turn the meter on by pressing the ON key.
- 4.Immerse the ph meter electrode in the solution to be tested.
- 5.Stir gently and wait for the reading to stabilize.
- 6.After use, rinse the electrode with water to minimize contamination. Turn the meter on by pressing the OFF key.
- 7.Always replace the protective cap after use.

lab instruments

PH Calibration

- 1.Remove the protective cap of electrode.
- 2.First rinse the electrode with distilled water, and suck it with filter-paper.
- 3.Turn the meter on by pressing the ON key.
- 4.Immerse the ph meter electrode in the solution to be tested.
- 5.Stir gently and wait for the reading to stabilize.
- 6.After use, rinse the electrode with water to minimize contamination. Turn the meter on by pressing the OFF key.

7. Always replace the protective cap after use.

lab instruments

Important:

The instrument's pH range must be re-calibrated whenever:

The electrode has been used (laid) for long time from the last calibration.

The electrode has been used in particularly taxing conditions.

The utmost accuracy is required.

Replace the new electrode.

lab instruments