

FGM-LX Binocular Gem Microscope /Gem Stereo Microscope /Stereo Zoom Microscope

Binocular stereo zoom microscope

The scope of application:

It can be used for identifying jewels and other precious stones.

Specifications:

Vertical Stage Microscope

Stereo zoom optics

180 degrees revolving objective

Zoom range:10-30X .Working distance:55mm

Superbright white LED source

Three kinds of illumination

Work normally in 100V-240V

Darkfield and brightfield controlled by electronic switching,easy to operate.

XYZ three-dimensional Motion Platform